

Happy

in the knowledge
that it's a
Swann Morton

When quality and precision count Dental Professionals count on Swann Morton.

With dental and facial aesthetics in mind Swann Morton offer a versatile range of handles and blades to meet the most demanding and complex requirements presented by modern day Periodontal, Cosmetic, Oral Implantology and Maxillofacial procedures.

When quality and precision count you can rely on Swann Morton for a durable blade that provides a consistent feel and sharpness resulting in better depth control from the first cut to the last.

Swann-Morton[®]
SHEFFIELD ENGLAND



**UK
CA**

All medical devices comply with the Medical Device Directive, 93/42/EEC as referred by Part II of the Medical Devices Regulations, 2002. Annex II excluding section 4 (as modified by Part 2 of Schedule 24 to the Medical Devices Regulations 2002).



EMS 633531



ENMS 633525



FM 73368




Scan this code for more product information

Swann-Morton®

ONE OF THE WORLD'S LEADING SURGICAL BLADE MANUFACTURERS


In addition to the more familiar traditional products Swann-Morton also provide a portfolio of blade and scalpel safety solutions available alongside their training programme on best practice when handling surgical blades which can also be accessed via their website in an online version for distance learning.

The Sharps Directive compliant sterile retractable safety scalpel supports existing safety protocols whilst encouraging safe handling and disposal of sharps within the dental practice. Originally available in blade shapes, 10, 11P and 15 which are colour-coded for ease of recognition the range has now been extended to include the 12 and 15c. Unit boxes of 10 have also been launched in addition to the original boxes of 25.



The single use sterile blade remover assists in the safe removal of a contaminated surgical blade from the traditional metal surgical handles as does the Clicksmart "BladeFLASK EVO" which can be conveniently mounted on each surgery wall and replaced when full.

The "Cygnetic" range combines a thicker stainless steel blade with a special innovative handle which grips the blade firmly on both sides providing additional lateral strength. The system minimises the handling of the blade when fitting and removing from the handle.



Let the quality, consistency and durability of Swann-Morton surgical blades and scalpels assist your practice in maximising efficiencies whilst delivering safe and successful outcomes for your patients.

When only Swann-Morton will do

Swann-Morton surgical blades can be used in the surgery for the following:

- Minor surgery
- Oral and Maxillofacial surgery
- Implantology
- Labial fraenum
- Incising abscesses
- Periodontal Surgery
- Maxillary sinus floor augmentation
- Raising Skin Flaps

More information can be found on the Swann-Morton website at www.swann-morton.com/dentist.php

Swann-Morton and the 'ring pattern' logo are registered trademarks of Swann-Morton Limited and related companies. The products shown are manufactured for use by trained health care professionals and can cause injury if used or handled incorrectly. The information and illustrations used in connection with the products shown on this leaflet are for general guidance only and users must seek proper training and instruction in their use before using or handling such products. Swann-Morton Limited exclude all liability for injury to the extent permitted by law.